

Product Catalogue
For
Food Hygiene and
Environmental Testing

HyServe GmbH & Co. KG

Table of Contents

1. Food Hygienic Test (Stamp Method):

<u>ID-No</u>	<u>Product Name</u>	<u>Application</u>	<u>page</u>
1000029	Food Stamp Standard Method Agar (100plates)	Total viable count	4
1000030	Food Stamp Standard Method Agar (30 plates)	Total viable count	4
1000041	Food Stamp X-GAL Agar (100 plates)	Coliforms	5
1000042	Food Stamp X-GAL Agar (30 plates)	Coliforms	5
1000043	Food Stamp XM-G Agar (100 plates)	<i>E. coli</i> / Coliforms	6
1000044	Food Stamp XM-G Agar (30 plates)	<i>E. coli</i> / Coliforms	6
1000031	Food Stamp TCBS Agar (100 plates)	<i>Vibrio parahaemolyticus</i>	7
1000032	Food Stamp TCBS Agar (30 plates)	<i>Vibrio parahaemolyticus</i>	7
1000033	Food Stamp TGSE Agar (100 plates)	<i>Staphylococcus aureus</i>	8
1000034	Food Stamp TGSE Agar (30 plates)	<i>Staphylococcus aureus</i>	8
1000039	Food Stamp MLCB Agar (100 plates)	<i>Salmonella</i>	9
1000040	Food Stamp MLCB Agar (30 plates)	<i>Salmonella</i>	9
1000035	Food Stamp Cereus Agar (100 plates)	<i>Bacillus cereus</i>	10
1000036	Food Stamp Cereus Agar (30 plates)	<i>Bacillus cereus</i>	10
1000027	Food Stamp Sabouraud Agar (100 plates)	Fungi	11
1000028	Food Stamp Sabouraud Agar (30 plates)	Fungi	11
1000037	Food Stamp Potato dextrose Agar with CP (100plates)	Food poisoning Fungi	12
1000038	Food Stamp Potato dextrose Agar with CP (30plates)	Food poisoning Fungi	12
1000046	Food Stamp XSA Agar (100 plates)	<i>Staphylococcus aureus</i>	13
1000047	Food Stamp XSA Agar (30 plates)	<i>Staphylococcus aureus</i>	13
1000025	Food Stamp XSA Agar (1000 plates)	<i>Coliforme</i>	14
1000026	Food Stamp XSA Agar (30 plates)	<i>Coliforme</i>	14

2. Environmental Test (Stamp Method)

<u>ID-No</u>	<u>Product Name</u>	<u>Application</u>	<u>page</u>
<u>For Medical Environment</u>			
06780	Clean Stamp MSO Agar (100 plates)	For MRSA	15
06781	Clean Stamp MSO Agar (30 plates)	For MRSA	15
1000179	Clean Stamp MSEY Agar (100 plates)	For Staphylococcus	16
1000180	Clean Stamp MSEY Agar (30 plates)	For Staphylococcus	16
<u>For Regular Environment</u>			
1000181	Clean Stamp SCD Agar (100 plates)	For Total viable count	17
1000182	Clean Stamp SCD Agar (30 plates)	For Total viable count	17
1000183	Clean Stamp SCDLP Agar (100 plates)	For Total viable count containing inactivator	18
1000184	Clean Stamp SCDLP Agar (30 plates)	For Total viable count containing inactivator	18
1000195	Clean Stamp25 SCD Agar (150 plates)	For Total viable count	19
1000196	Clean Stamp25 SCD Agar (30 plates)	For Total viable count	19
1000197	Clean Stamp25 SCDLP Agar (150 plates)	For Total viable count containing inactivator	20
1000198	Clean Stamp25 SCDLP Agar (30 plates)	For Total viable count containing inactivator	20
1000199	Clean Stamp25 CPSB Agar (150 plates)	For Fungi	21
1000200	Clean Stamp25 CPSB Agar (30 plates)	For Fungi	21

3. Food Allergen Detection Kit

<u>ID-No</u>	<u>Product Name</u>	<u>Application</u>	<u>page</u>
08600	FASTKIT ELISA Egg Kit	For Food Allergen screening (Egg)	22
08601	FASTKIT ELISA Milk Kit	For Food Allergen screening (Milk)	22
08602	FASTKIT ELISA Wheat Kit	For Food Allergen screening (Wheat)	22
08603	FASTKIT ELISA Buckwheat	For Food Allergen screening (Buckwheat)	22
08604	FASTKIT ELISA Peanut	For Food Allergen screening (Peanut)	22

4. Culture Media (dehydrated)

ID-No.	Product Name	Application	page
<u>For Salmonella</u>			
1000071	EEM Broth (100g)	Enrichment broth for <i>Salmonella</i>	23
05131	Buffered Peptone Water (BPW) (300g)	Enrichment broth for <i>Salmonella</i>	24
05021	SS Agar (4L)	Isolation of <i>Salmonella</i>	25
05020	SS Agar (20L)	Isolation of <i>Salmonella</i>	25
05032	SS Agar with Sucrose (280g)	Isolation of <i>Salmonella</i>	26
05033	SS Agar with Sucrose (20L)	Isolation of <i>Salmonella</i>	26
<u>For E. coli / Coliforms</u>			
1000211	Blue Light Broth (300g)	Rapid detection of <i>E. coli</i> and Coliforms from Water and Foods	27
05036	MacConkey Agar (300g)	Isolation of Enteric bacteria	28
05643	MacConkey Sorbitol Agar (300g)	Isolation for <i>E. coli</i> O157	29
05639	Lauryl Sulfate MUG Broth (300g)	Rapid detection of <i>E. coli</i> and Coliforms from Water and Foods	30
05642	X-GAL Agar (300g)	Chromogenic enzyme substrate agar for Coliform isolation	31
05647	XM-G Agar (300g)	Chromogenic enzyme substrate agar for <i>E. coli</i> and Coliforms	32
<u>For Staphylococcus</u>			
05234	Staphylococcus Medium No. 110 (300g)	Isolation of <i>Staphylococcus</i>	33
05236	Mannitol Salt Agar (300g)	Isolation of <i>Staphylococcus</i>	34
<u>For Enterococci</u>			
05679	EF Agar Base (100g)	Isolation of Enterococci from Water, Milk and Meat products	35
<u>For Bacillus cereus</u>			
05282	NGKG Agar Base (300g)	Selective isolation of <i>B. cereus</i> from Foods	36
<u>For Anaerobes</u>			
05426	GAM Agar, Modified (300g)	Isolation for Anaerobic bacteria	37
05433	GAM Broth, Modified (100g)	Isolation for Anaerobic bacteria	38
05440	<i>Bacteroides</i> Agar (100g)	Differentiation and isolation for <i>Bacteroides</i>	39
05441	FM Agar (100g)	Differentiation and isolation for Fusobacterium	40
06593	Anaero Mate-P (50 sets for 100 plates)	Simple bag for Anaerobic culture	41
06594	Anaero Mate-J (50 sets for 500 plates)	Simple bag for Anaerobic culture	41
<u>For general use</u>			
05503	Heart Infusion Agar (300g)	Isolation for various bacteria	42
05505	Heart Infusion Broth (100g)	Isolation for various bacteria	43
05506	Brain Heart Infusion Agar (300g)	Isolation for various bacteria	44
05508	Brain Heart Infusion Broth (300g)	Isolation for various bacteria	45
05516	Trypto-Soya Agar (SCD Agar) (300g)	Isolation for various bacteria	46
05533	Mueller-Hinton Agar-N (300g)	Isolation for various bacteria	47
<u>For Yeast and Mold</u>			
05701	Sabouraud Agar (300g)	Isolation of Yeast and Mold	48
05703	Candida GE Agar (100g)	Isolation for Candida	49

5. Identification Test Kit

ID-No.	Product Name	Application	page
06626	ID Test EB-20 (25 tests)	Identification test for Enterobacteriaceae	50
06629	ID Test NF-18 (25 tests)	Identification test for Gram Negative Non-fermenters	51
06637	ID Test SP-18 (25 tests)	Identification test for Staphylococcus	52

Food Stamp Standard Method Agar (SMA)

Simple and Easy Stamp Medium for Food Hygiene test: Total Viable Bacterial Count

Code **1000029 – 100 plates**
 1000030 – 30 plates

Formula (in 1 liter)

Peptone.....	5.0g	Yeast Extract.....	2.5g
Dextrose.....	1.0g	Agar.....	15.0g

Directions

Food Stamp is a prepared agar medium for Stamp method, on which agar stands up slightly above the rim of special Petri dish of 10 cm². Take off the cap of Food Stamp and gently press the medium against the surface of specimen. The surface of agar is elastic enough to be pressed firmly against the specimen. Press against the different parts of the specimen when several kinds of Food Stamps are tested simultaneously. Put the cap again immediately after pressing.

This Food Stamp is designed to measure and detect the degree of contamination of the specimen.

Incubate at 37°C for 24 - 48 hours or at Room Temperature for 48 – 96 hours.

Interpretations

Count all colonies grown on the surface.

Storage

Keep at 4 – 10°C. Do not freeze.
Six (6) months after manufacturing.

Further Information

See Page 53



Food Stamp X-GAL Agar

Simple and Easy Stamp Medium for Food Hygiene test: Coliforms

Code 1000041 – 100 plates
1000042 – 30 plates

Formula (in 1 liter)

Peptone.....	15.0g	Yeast Extract.....	5.0g
Sodium Pyruvate.....	1.0g	Sodium Chloride.....	15.0g
Disodium Phosphate.....	2.0g	Potassium Nitrate.....	1.0g
Sodium Lauryl Sulfate.....	0.15g	Agar.....	15.0g
5-bromo-4-chloro-3-indolyl- β -D-galactopyranoside (X-GAL).....		0.15g	

Directions

Food Stamp is a prepared agar medium for Stamp method, on which agar stands up slightly above the rim of special Petri dish of 10 cm². Take off the cap of Food Stamp and gently press the medium against the surface of specimen. The surface of agar is elastic enough to be pressed firmly against the specimen. Press against the different parts of the specimen when several kinds of Food Stamps are tested simultaneously. Put the cap again immediately after pressing. Incubate at 37°C for 24 - 48 hours.

Coliform group decompose X-GAL (colorimetric enzyme substrate) in the medium to bring out blue / blue-green color by β -galactosidase produced by Coliforms.

Growth of all other bacteria will be inhibited, and they develop only white colonies even are they grow.

Interpretations

Count all blue / blue-green color colonies grown on the surface.


Storage

Keep at 4 – 10°C. Do not freeze.

Five (5) months after manufacturing.

Further Information

See Page 53



Food Stamp XM-G Agar

Simple and Easy Stamp Medium for Food Hygiene test: *Escherichia coli* and Coliforms

Code 1000043 – 100 plates
1000044 – 30 plates

Formula (in 1 liter)

Peptone.....	15.0g	L-Tryptophan.....	1.0g
Sodium Pyruvate.....	1.0g	D-Sorbitol.....	1.0g
Sodium Chloride.....	5.0g	Monopotassium Phosphate.....	2.2g
Dipotassium Phosphate.....	2.7g	Potassium Nitrate.....	1.0g
Sodium Lauryl Sulfate.....	0.2g	Agar.....	15.0g
5-bromo-4-chloro-3-indolyl- β -D-glucuronide (X-GLUC).....	0.1g		
5-bromo-6-chloro-3-indolyl- β -D-galactopyranoside (MAGENTA-GAL)	0.1g		

Directions

Food Stamp is a prepared agar medium for Stamp method, on which agar stands up slightly above the rim of special Petri dish of 10 cm². Take off the cap of Food Stamp and gently press the medium against the surface of specimen. The surface of agar is elastic enough to be pressed firmly against the specimen. Press against the different parts of the specimen when several kinds of Food Stamps are tested simultaneously. Put the cap again immediately after pressing. Incubate at 37°C for 24 - 48 hours. Coliform group decompose X-GAL (colorimetric enzyme substrate) in the medium to bring out blue / blue-green color by β -galactosidase produced by Coliforms. Growth of all other bacteria will be inhibited, and they develop only white colonies even are they grow.

Interpretations

Count all pink / red-purple color colonies grown on the surface as Coliforms.

Count all blue / blue-purple color colonies grown on the surface as *E. coli*.

E. coli O-157 does not have a β -glucuronidase, and then it will be identified as a Coliform.

Overtime incubation may foster growth of microorganisms other than *E. coli* and Coliforms.

Red color may be observed if the sample contains lactobacilli that has also a β -galactosidase.


Storage

Keep at 4 – 10°C. Do not freeze.

Five (5) months after manufacturing.

Further Information

See Page 53



Food Stamp TCBS Agar

Simple and Easy Stamp Medium for Food Hygiene test: *Vibrio parahaemolyticus*

Code 1000031 – 100 plates
1000032 – 30 plates

Formula (in 1 liter)

Peptone.....	10.0g	Ferric Citrate.....	1.0g
Yeast Extract.....	1.0g	Saccharose.....	17.0g
Sodium Chloride.....	10.0g	Ox gall.....	5.0g
Sodium Thiosulfate.....	10.0g	Bromthymol Blue.....	0.04g
Sodium Citrate.....	10.0g	Thymol Blue.....	0.04g
Sodium Cholate.....	3.0g	Agar.....	15.0g

Directions

Food Stamp is a prepared agar medium for Stamp method, on which agar stands up slightly above the rim of special Petri dish of 10 cm². Take off the cap of Food Stamp and gently press the medium against the surface of specimen. The surface of agar is elastic enough to be pressed firmly against the specimen. Press against the different parts of the specimen when several kinds of Food Stamps are tested simultaneously. Put the cap again immediately after pressing.

Incubate at 37°C for 24 - 48 hours.

Interpretations


Green (*V. parahaemolyticus*) and Yellow (*V. alginolyticus*) colonies are observed on the surface. Possible contamination of *V. parahaemolyticus* is suspected when many yellow colonies (*V. alginolyticus*) are observed.

Storage

Keep at 4 – 10°C. Do not freeze.
Six (6) months after manufacturing.

Further Information

See Page 53



Food Stamp TGSE Agar

Simple and Easy Stamp Medium for Food Hygiene test: *Staphylococcus aureus*

Code 1000033 – 100 plates
1000034 – 30 plates

Formula (in 1 liter)

Peptone.....	15.0g	Soya Peptone.....	2.0g
Beef Extract.....	3.0g	Calcium Chloride.....	17.0g
Sodium Chloride.....	65.0g	Mannitol.....	5.0g
Lithium Chloride . .	10.0g	Glycine.....	0.04g
Potassium Tellurite.....	10.0g	Egg yolk.....	5%
Agar.....	15.0g		

Directions

Food Stamp is a prepared agar medium for Stamp method, on which agar stands up slightly above the rim of special Petri dish of 10 cm². Take off the cap of Food Stamp and gently press the medium against the surface of specimen. The surface of agar is elastic enough to be pressed firmly against the specimen. Press against the different parts of the specimen when several kinds of Food Stamps are tested simultaneously. Put the cap again immediately after pressing.

Incubate at 37°C for 24 - 48 hours or at Room Temperature for 72 – 96 hours.

Interpretations


Staphylococcus aureus forms black colonies with milky surroundings around the colonies. *S. aureus* shows a positive egg yolk reaction. Black colony with negative egg yolk reaction is not interpreted as *S. aureus*.

Storage

Keep at 4 – 10°C. Do not freeze.
Six (6) months after manufacturing.

Further Information

See Page 53



Food Stamp MLCB Agar

Simple and Easy Stamp Medium for Food Hygiene test: *Salmonella*

Code 1000039 – 100 plates
1000040 – 30 plates

Formula (in 1 liter)

Peptone.....	10.0g	Sodium Thiosulfate.....	4.0g
Yeast Extract.....	3.0g	Ferric Ammonium Citrate.....	1.0g
Heart Extract Powder.	2.0g	Brilliant Green.....	0.0125g
Sodium Chloride.	4.0g	Crystal Violet.....	0.01g
Mannitol.....	3.0g	Agar.....	15.0g
L-Lysine Hydrochloride.....	5.0g		

Directions

Food Stamp is a prepared agar medium for Stamp method, on which agar stands up slightly above the rim of special Petri dish of 10 cm². Take off the cap of Food Stamp and gently press the medium against the surface of specimen. The surface of agar is elastic enough to be pressed firmly against the specimen. Press against the different parts of the specimen when several kinds of Food Stamps are tested simultaneously. Put the cap again immediately after pressing.

Incubate at 37°C for 24 - 48 hours.

Interpretations

Salmonella produces hydrogen sulfide and forms black colonies or colonies with black center.

Citrobacter may develop black colonies just like *Salmonella*.

Purple colonies are not *Salmonella*.


Storage

Keep at 4 – 10°C. Do not freeze.

Four (4) months after manufacturing.

Further Information

See Page 53



Food Stamp Cereus Agar

Simple and Easy Stamp Medium for Food Hygiene test: *Bacillus cereus*

Code 1000035 – 100 plates
1000036 – 30 plates

Formula (in 1 liter)

Proteose Peptone.....	10.0g	Glycine.....	10.0g
Peptone.....	10.0g	Phenol Red.....	0.05g
Heart Extract Powder.	5.0g	Polymyxin B.....	50,000unit
Sodium Chloride	5.0g	Egg Yolk Suspension.....	10%
Lactose	10.0g	Agar.....	25.0g

Directions

Food Stamp is a prepared agar medium for Stamp method, on which agar stands up slightly above the rim of special Petri dish of 10 cm². Take off the cap of Food Stamp and gently press the medium against the surface of specimen. The surface of agar is elastic enough to be pressed firmly against the specimen. Press against the different parts of the specimen when several kinds of Food Stamps are tested simultaneously. Put the cap again immediately after pressing.

Incubate at 37°C for 24 - 48 hours.

Interpretations

Bacillus cereus forms white colonies with an irregular rim, which develop opaque zone (positive egg yolk reaction) around the colonies and change the color of medium to Red.
Small colonies with negative egg yolk reaction are not *B. cereus*.

Storage

Keep at 4 – 10°C. Do not freeze.
Six (6) months after manufacturing.

Further Information

See Page 53



Food Stamp Sabouraud Agar

Simple and Easy Stamp Medium for Food Hygiene test: *Bacillus cereus*

Code 1000027 – 100 plates
1000028 – 30 plates

Formula (in 1 liter)

Peptone.....	10.0g	Dextrose.....	40.0g
Agar.....	15.0g		

Directions

Food Stamp is a prepared agar medium for Stamp method, on which agar stands up slightly above the rim of special Petri dish of 10 cm². Take off the cap of Food Stamp and gently press the medium against the surface of specimen. The surface of agar is elastic enough to be pressed firmly against the specimen. Press against the different parts of the specimen when several kinds of Food Stamps are tested simultaneously. Put the cap again immediately after pressing.

Incubate at 30°C for 48 - 72 hours.

Interpretations

Fungi develop characteristic fluffy colonies on the surface.
All characteristic colonies should be counted for evaluation.

Storage

Keep at 4 – 10°C. Do not freeze.
Six (6) months after manufacturing.

Further Information

See Page 53



Food Stamp Potato Dextrose Agar with CP

Simple and Easy Stamp Medium for Food Hygiene test: Food poisoning Fungi

Code 1000037 – 100 plates
1000038 – 30 plates

Formula (in 1 liter)

Potato Extract.....	4.0g	Dextrose.....	20.0g
Chloramphenicol.....	100mg	Agar.....	15.0g

Directions

Food Stamp is a prepared agar medium for Stamp method, on which agar stands up slightly above the rim of special Petri dish of 10 cm². Take off the cap of Food Stamp and gently press the medium against the surface of specimen. The surface of agar is elastic enough to be pressed firmly against the specimen. Press against the different parts of the specimen when several kinds of Food Stamps are tested simultaneously. Put the cap again immediately after pressing.

Incubate at 30°C for 48 - 72 hours.

Interpretations


Since Chloramphenicol in the medium inhibits the growths of other bacteria, all colonies can be counted for Fungi.

Storage

Keep at 4 – 10°C. Do not freeze.
Five (5) months after manufacturing.

Further Information

See Page 53



Food Stamp X-SA Agar

Simple and Easy Stamp Medium for Food Hygiene test: *Staphylococcus aureus*

Code 1000046 – 100 plates
1000047 – 30 plates

Formula (in 1 liter)

Peptone	13.0g	Beef Extract.....	3.0g
Lithium Chloride.....	5.0g	Mannitol.....	10.0g
Agar.....	14.0g	Chloramphenicol.....	100mg
5-bromo-4-chloro-3-indolyl-phosphate,disodium salt sesquihydrate (X-SA).....			

Directions

Food Stamp is a prepared agar medium for Stamp method, on which agar stands up slightly above the rim of special Petri dish of 10 cm². Take off the cap of Food Stamp and gently press the medium against the surface of specimen. The surface of agar is elastic enough to be pressed firmly against the specimen. Press against the different parts of the specimen when several kinds of Food Stamps are tested simultaneously. Put the cap again immediately after pressing.

Incubate at 35 - 37°C for 22 - 24 hours.

Interpretations

Staphylococcus aureus forms blue color smooth colonies with converging lens shape colonies.

*Coagulase-negative *Staphylococcus* form white or blue smallness colonies.

**Bacillus* sp. may form light blue and flat colonies without glaze on this plate.


Storage

Keep at 4 – 10°C. Do not freeze.

Four (4) months after manufacturing.

Further Information

See Page 53



Food Stamp Standard Desoxycholat Agar (DESO)

Simple and Easy Stamp Medium for Food Hygiene test: Coliforme

Code 1000025 – 100 plates
1000026 – 30 plates

Formula (in 1 liter)

Sodium Desoxycholate	1.0g	Dipotassium Phosphate	2.0g
Sodium Chlorid	5.0g	Ferric Ammonium Citrate	2.0g
Peptone	10.0g	Lactose	10.0g
Neutral Red	0.033g	Agar	15.0g

Directions

Food Stamp is a prepared agar medium for Stamp method, on which agar stands up slightly above the rim of special Petri dish of 10 cm². Take off the cap of Food Stamp and gently press the medium against the surface of specimen. The surface of agar is elastic enough to be pressed firmly against the specimen. Press against the different parts of the specimen when several kinds of Food Stamps are tested simultaneously. Put the cap again immediately after pressing.

Incubate at 37°C for 24 - 48 hours or room temperature for 48 hours

Interpretations


Count all colonies on the surface

Storage

Keep at 4 – 10°C. Do not freeze.
Four (4) months after manufacturing.

Further Information

See Page 53



Clean Stamp MSO Agar (MSO)

Environmental Test for Medical Ambience: MRSA

Code 06780 – 100 plates
06781 – 30 plates

Formula (in 1 liter)

Beef Extract.....	1.0g	Colistin.....	10.0mg
Peptone.....	10.0g	Amphotericin B.....	2.0mg
Sodium Chloride.....	40.0g	Aztreonam.....	5.0mg
Mannitol.....	10.0g	Phenol Red.....	25.0mg
Sodium Pyruvate.....	2.0g	Egg Yolk.....	30.0ml
Lithium Chloride.....	5.0g	Agar.....	15.0g
Oxacillin.....	6.0mg		
pH 7.4 ± 0.1			

Directions

Clean Stamp is a prepared agar medium for Stamp method, on which agar stands up slightly above the rim of special Petri dish of 10 cm². Take off the cap of Clean Stamp and gently press the medium against the surface of specimen. The surface of agar is elastic enough to be pressed firmly against the specimen. Press against the different parts of the specimen when several kinds of Clean Stamps are tested simultaneously. Put the cap again immediately after pressing. Incubate at 35°C for 24 - 48 hours.

Interpretations

MRSA: MRSA forms round yellow/cream color colonies with yellow surroundings around the colonies. Due to positive Egg Yolk reaction, white turbidity is developed inside the medium around the colonies while the surface has pearly luster.

In case of faint Egg Yolk reaction is observed, extend incubation to 48 hours, and read the reaction after colonies are removed.

Non-MRSA: Mannitol non-fermenting (no yellow color change) and negative Egg Yolk reaction bacteria are not MRSA.

Storage

Do keep at 4 – 10°C. Do not freeze.
Two (2) months after manufacturing.

Further Information

See Page 53

Clean Stamp Mannitol Salt Agar with Egg Yolk (MSEY)

Microbiological Test for Medical Ambience: *Staphylococcus*

Code 1000179 – 100 plates
 1000180 – 30 plates

Formula (in 1 liter)

Beef Extract.....	1.0g	Phenol Red.....	25.0mg
Peptone.....	10.0g	Egg Yolk.....	30.0ml
Sodium Chloride.....	75.0g	Agar.....	15.0g
Mannitol.....	10.0g		
pH 7.4 ± 0.1			

Directions

Clean Stamp is a prepared agar medium for Stamp method, on which agar stands up slightly above the rim of special Petri dish of 10 cm². Take off the cap of Clean Stamp and gently press the medium against the surface of specimen. The surface of agar is elastic enough to be pressed firmly against the specimen. Press against the different parts of the specimen when several kinds of Clean Stamps are tested simultaneously. Put the cap again immediately after pressing. Incubate at 35°C for 24 - 48 hours.

Interpretations

***Staphylococcus aureus*:** *Staphylococcus aureus* forms round yellow/cream color colonies with yellow surroundings around the colonies. Due to positive Egg Yolk reaction, white turbidity is developed inside the medium around the colonies while the surface has pearly luster. In case of faint Egg Yolk reaction is observed, extend incubation to 48 hours, and read the reaction after colonies are removed.

Other *Staphylococcus*: Other *Staphylococcus* forms small round white colonies with no color changes or red surroundings around the colonies, and do not show any positive Egg Yolk reaction.

Storage

Keep at 4 – 10°C. Do not freeze.
Five (5) months after manufacturing.

Further Information

See Page 53

Clean Stamp SCD Agar (SCD)

Microbiological Test for Environment: Viable Bacterial Count
Soybean Casein Digest Agar

Code 1000181 – 100 plates
1000182 – 30 plates

Formula (in 1 liter)

Peptone.....	15.0g	Sodium Chloride.....	5.0g
Soya Peptone.....	5.0g	Agar.....	15.0g
pH 7.2±0.1			

Directions

Clean Stamp is a prepared agar medium for Stamp method, on which agar stands up slightly above the rim of special Petri dish of 10 cm². Take off the cap of Clean Stamp and gently press the medium against the surface of specimen. The surface of agar is elastic enough to be pressed firmly against the specimen. Press against the different parts of the specimen when several kinds of Clean Stamps are tested simultaneously. Put the cap again immediately after pressing. Incubate at 35°C for 24 - 48 hours.

Interpretations

Count all colonies grown on the surface.

Storage

Keep at 4 – 10°C. Do not freeze.
Five (5) months after manufacturing.

Further Information

See Page 53

Clean Stamp

SCD Agar with Lecithin, Polysorbate 80 (SCDLP)

Microbiological Test for Environment: Viable Bacterial Count
Including inactivate agents

Code **1000183 – 100 plates**
 1000184 – 30 plates

Formula (in 1 liter)

Peptone.....15.0g	Lecithin.....1.0g
Soya Peptone.....5.0g	Polysorbate 80.....7.0g
Sodium Chloride.....5.0g	Agar.....15.0g
pH 7.2±0.1	

Directions

Clean Stamp is a prepared agar medium for Stamp method, on which agar stands up slightly above the rim of special Petri dish of 10 cm². Take off the cap of Clean Stamp and gently press the medium against the surface of specimen. The surface of agar is elastic enough to be pressed firmly against the specimen. Press against the different parts of the specimen when several kinds of Clean Stamps are tested simultaneously. Put the cap again immediately after pressing.

This Clean Stamp is designed to measure and detect the degree of pollution of the specimen which surface has been processed with chemicals or disinfectants. As lecithin and polysorbate 80 inactivate those chemicals and disinfectants absorbed into medium, more reliable results can be obtained.

Incubate at 35°C for 24 - 48 hours.

Interpretations

Count all colonies grown on the surface.

Storage

Keep at 4 – 10°C. Do not freeze.
Five (5) months after manufacturing.

Further Information

See Page 53

Clean Stamp 25 SCD Agar (SCD)

Microbiological Test for Environment: Viable Bacterial Count
Soybean Casein Digest Agar

Code 1000196 – 150 plates
1000195 – 30 plates

Formula (in 1 liter)

Peptone.....	15.0g	Sodium Chloride.....	5.0g
Soya Peptone.....	5.0g	Agar.....	15.0g
pH 7.2±0.1			

Directions

Clean Stamp 25 is a prepared agar medium for Stamp method, on which agar stands up slightly above the rim of special Petri dish of 25 cm². Take off the cap of Clean Stamp 25 and gently press the medium against the surface of specimen. The surface of agar is elastic enough to be pressed firmly against the specimen. Press against the different parts of the specimen when several kinds of Clean Stamps are tested simultaneously. Put the cap again immediately after pressing.

Incubate at 30 - 35°C for 2 – 5 days.

Interpretations

Count all colonies grown on the surface.

Remarks

Clean Stamp 25 SCD is made and prepared of Code 05516 Trypto-Soya Agar according to its directions, and aseptically distributed and solidified in the special Petri dishes having an inside diameter of 58mm, outside diameter of 73mm and height of 18mm. One clean film bag contains 3 Clean Stamp 25s connected in series.

Storage

Keep at 4 – 10°C. Do not freeze.

Five (5) months after manufacturing.

Further Information

See Page 53

Clean Stamp 25 SCDLP Agar (SCDLP)

Microbiological Test for Environment: Viable Bacterial Count
Including inactivate agents

Code **1000197 – 150 plates**
 1000198 – 30 plates

Formula (in 1 liter)

Peptone.....15.0g	Lecithin.....1.0g
Soya Peptone.....5.0g	Polysorbate 80.....7.0g
Sodium Chloride.....5.0g	Agar.....15.0g
pH 7.2±0.1	

Directions

Clean Stamp 25 SCDLP is a prepared agar medium for Stamp method, on which agar stands up slightly above the rim of special Petri dish of 25 cm². Take off the cap of Clean Stamp 25 SCDLP and gently press the medium against the surface of specimen. The surface of agar is elastic enough to be pressed firmly against the specimen. Press against the different parts of the specimen when several kinds of Clean Stamps are tested simultaneously. Put the cap again immediately after pressing.

This Clean Stamp 25 SCDLP is designed to measure and detect the degree of pollution of the specimen which surface has been processed with chemicals or disinfectants. As lecithin and polysorbate 80 inactivate those chemicals and disinfectants absorbed into medium, more reliable results can be obtained.

Incubate at 30 - 35°C for 2 – 5 days.

Interpretations

Count all colonies grown on the surface.

Remarks

Clean Stamp 25 SCDLP is made and prepared of Code 05516 Trypto-Soya Agar Plus Lecithin and Polysorbate 80 according to its directions, and aseptically distributed and solidified in the special Petri dishes having an inside diameter of 58mm, outside diameter of 73mm and height of 18mm. One clean film bag contains 3 Clean Stamp 25s connected in series.

Storage

Keep at 4 – 10°C. Do not freeze.

Five (5) months after manufacturing.

Further Information

See Page 53

Clean Stamp 25
Sabouraud Agar with Chloramphenicol
(CPSB)

Microbiological Test for Environment: Fungus

Code **1000199 – 150 plates**
 1000200 – 30 plates

Formula (in 1 liter)

Peptone.....	10.0g	Chloramphenicol.....	50mg
Dextrose.....	40.0g	Agar.....	15.0g
pH 5.9±0.2			

Directions

Clean Stamp 25 CPSB is a prepared agar medium for Stamp method, on which agar stands up slightly above the rim of special Petri dish of 25 cm². Take off the cap of Clean Stamp 25 CPSB and gently press the medium against the surface of specimen. The surface of agar is elastic enough to be pressed firmly against the specimen. Press against the different parts of the specimen when several kinds of Clean Stamps are tested simultaneously. Put the cap again immediately after pressing.

Interpretations

Count all colonies grown on the surface.

Remarks

Clean Stamp 25 CPSB is made and prepared of Code 05701 Sabouraud Agar Plus Chloramphenicol according to its directions, and aseptically distributed and solidified in the special Petri dishes having an inside diameter of 58mm, outside diameter of 73mm and height of 18mm. One clean film bag contains 3 Clean Stamp 25s connected in series.

Storage

Keep at 4 – 10°C. Do not freeze.
Five (5) months after manufacturing.

Further Information

See Page 53

Food Allergen Screening Test Kit FASTKIT Series

For raw materials like Egg, Milk, Wheat, Buckwheat (Soba) and Peanut

Code 08600 - FASTKIT ELISA Egg
08601 - FASTKIT ELISA Milk
08602 - FASTKIT ELISA Wheat
08603 - FASTKIT ELISA Buckwheat (Soba)
08604 - FASTKIT ELISA Peanut

Directions

The regulation related Food Sanitation Law in Japan was revised, and new regulation of Labeling of Allergenic Substance in Food was issued to make it compulsory to indicate on the label of the name of specific raw materials, its interfusion or carry-over. Specific raw materials stipulated are egg, milk, wheat, soba and peanut, all of which may be considered to be allergenic foods and may cause serious allergy.

FAST (Food Allergen Screening Test) KIT series are the test to detect the proteins included in the specific raw materials, and designed for all those five raw materials in both rude and processed materials specified by regulation.

Features and Benefit

1. **High sensitivity:** FASTKIT can detect and measure as low as $1 \sim 5 \mu\text{g/g}$ (ppm) of specific raw material mixed in food.
It is said that allergic reactions may be frequently caused in case of the foods contained more than several dozen $\mu\text{g/g}$ of antigenic proteins in them. FASTKIT can detect and measure a small amount of specific raw material mixed in food ($1 \sim 5 \mu\text{g/g}$) depending on the foods.
2. **Polyclonal antibody that can detect several antigenic proteins in the food all together.**
Unlike conventional methods that detect only singly antigenic protein, FASTKIT consists of specific antibody that detects several antigenic proteins in specific raw materials. For example, the test consisted of the antibody that reacts only with egg white (ex. OVA), cannot detect egg yolk in the food, and no indication is eventually made for the egg yolk, which violate the regulation of Food Sanitation Law. Instead FASTKIT detects both of them.
3. **Well suited even for processed foods.**
FASTKIT consists of antibody that is immunized even with heat-denatured proteins, and FASTKIT can detect specific raw materials in not only rude but also foods processed through heat or pressure during manufacturing.

Measurement Principle

ELISA Method (enzyme-linked immunosorbent assay)

- Detectable level (Standard solution): $1\text{ng/mL} \sim 100\text{ng/mL}$
- Reproducibility (CV value): Less than 10%

Operating Procedure with food sample (example)

1. Food sample
 2. Add 19 fold of extraction buffer solution, and extract after crushing the sample with homogenizer
 3. Centrifuge, and remove insoluble matters through filtration
 4. Dilute to 10 fold with dilution buffer solution
 5. ELISA Testing (3.5 hours)
 6. Read the results with Plate reader
- Egg, Milk and Wheat Kit: 450nm, Buckwheat (Soba) and Peanut Kit: 405nm
 - Extraction method is subject to foodstuff.

Kit Composition

A: Microtiter plate.....	96 well x 1 plate
B: Standard solution.....	500 μL x 1 bottle
C: Buffer solution for diluting.....	100mL x 1 bottle
D: Antibody conjugated with biotin.....	150 μL x 1 bottle
E: Enzyme-avidin binding Solution.....	150 μL x 1 bottle
F: Color Forming reagent.....	12mL x 1 bottle
G: Concentrated Extraction Buffer.....	100mL x 1 bottle
H: Stopping Solution.....	12mL x 1 bottle
I: Concentrated Washing Solution.....	100mL x 1 bottle

Further Information

See Page 53

EEM Broth

Enterobacteriaceae Enrichment Mannitol Broth for
Pre-culture of *Salmonella*

Code 1000071 – 100g

Formula (in 43.5g for 1 liter)

Peptone.....	10.0g	Monopotassium Phosphate.....	2.0g
Mannitol.....	5.0g	Bile Powder.....	20.0g
Disodium Phosphate.....	6.45g	Brilliant Green.....	0.0135g
pH 7.2 ± 0.1			

Directions

Suspend 43.5g of the dehydrated medium in 1,000ml of distilled water, mix well to dissolve the medium. Distribute adequate amounts into sterilized test tubes or sterilized flask and then heat them at 100°C for 30 minutes. Avoid excessive heating, or do not autoclave. Add the specimen of about one-tenth of the medium and incubate at 37°C for 3 hours. Stir them well after the first incubation and continue the incubation for 20 – 24 hours or more. Finally, take a loopful of the culture and subculture on an adequate selective medium.

Remarks

The medium is quite suitable for preliminary culture of *Salmonella* in food and feed.

Storage

Keep dry at room temperature.
Three (3) years after manufacturing.

Further Information

See Page 53

Buffered Peptone Water (BPW)Buffered Peptone Water for Pre-culture of *Salmonella***Code 05131 – 300g****Formula (in 20.0g for 1 liter)**

Peptone.....	4.5g	Monopotassium Phosphate.....	1.5g
Sodium Chloride.....	7.2g	Disodium Phosphate.....	3.5g
pH 7.2±			

Directions

Suspend 20.0g of the dehydrated medium in 1,000ml of distilled water, mix well to dissolve the medium. Distribute adequate amounts into sterilized test tubes or sterilized flask and then autoclave at 121°C for 15 minutes.

Add 1/10 volume of the homogenized specimen (example: 1mL homogenized specimen put into 10mL broth) into the broth and stir them well. Incubate at 37°C for 20 – 24 hours. Finally, take one loop of the culture and subculture on an adequate selective medium (ex. Compact Dry SL, SS agar, MLCB Agar).

For liquid Egg sample, add 0.2g L-cystein or 64mg FeSO₄ • 7H₂O in 1 liter of Broth.

Remarks

The medium is suitable for pre-culture of *Salmonella* in food.

Storage

Keep dry at room temperature.

Three (3) years after manufacturing.

Further Information

See Page 53

SS Agar

Selective isolation and differential medium for
Salmonella and *Shigella*

Code 05021 – 4L (240g (granule))
05020 – 20L (1,200g (granule))

Formula (in 60.0g for 1 liter)

Beef Extract.....	5.0g	Sodium Thiosulfate.....	5.5g
Peptone.....	7.5g	Ferric Citrate.....	1.0g
Bile Salts.....	9.0g	Neutral Red.....	0.025g
Lactose.....	10.0g	Brilliant Green.....	0.00033g
Sodium Citrate.....	8.5g	Agar.....	13.5g
pH 7.3 ± 0.05			

Directions

Suspend 60.0g of the dehydrated medium in 1,000ml of distilled water, mix well and heat to dissolve the medium. Cool to about 50°C and distribute about 20ml amounts to Petri dishes. Dry the surface of the plate until moisture disappears from the surface, leaving the lids of the dishes partly open for about an hour.

The medium does not require autoclaving or any other aseptic procedures.

Determinations

Inoculate the specimens heavily and incubate at 37°C for 20 – 24 hours.

Pathogenic bacteria form semitransparent colonies

Lactose fermenters form pink or red colonies. Some strains of *Proteus* produce hydrogen sulfide and turn black. *Proteus* grown on the medium forms colonies similar to those of pathogenic bacteria.

Remarks

Features

1. The medium inhibits the growth of *E. coli* and miscellaneous bacteria.
2. The medium does not require any sterilizing procedures.
3. The medium enables to detect pathogenic bacteria easily and its detection rate is far more excellent than other media.

Note: The medium requires inoculation of a large amount of feces than in case of other selective media.

Storage

Keep dry at room temperature.

Three (3) years after manufacturing.

Further Information

See Page 53

SS Agar with Sucrose

Selective isolation and differential medium for
Salmonella and *Shigella*

Code 05032 – 280g (granule)
05033 – 20L (70g x 20 (granule))

Formula (in 70.0g for 1 liter)

Beef Extract.....	5.0g	Sodium Thiosulfate.....	5.5g
Peptone.....	8.0g	Ferric Citrate.....	1.0g
Bile Salts.....	8.5g	Neutral Red.....	0.025g
Lactose.....	10.0g	Brilliant Green.....	0.00033g
Sodium Citrate.....	8.5g	Agar.....	13.5g
Sucrose.....	10.0g		
pH 7.2± 0.05			

Directions

Suspend 60.0g of the dehydrated medium in 1,000ml of distilled water, mix well and heat to dissolve the medium. Cool to about 50°C and distribute about 20ml amounts to Petri dishes. Dry the surface of the plate until moisture disappears from the surface, leaving the lids of the dishes partly open for about one hour.

The medium does not require autoclaving or any other aseptic procedures.

Determinations

Inoculate the specimens heavily and incubate at 37°C for 20 – 24 hours.

Pathogenic bacteria form semitransparent colonies

Lactose fermenters form pink or red colonies. Also Sucrose fermented *Citrobacter* does not form black colonies and it can be identified easily.

Remarks

Features

- The medium inhibits the growth of *E. coli* and miscellaneous bacteria.
- The medium does not require any sterilizing procedures.
- The medium enables to detect pathogenic bacteria easily and its detection rate is far more excellent than other media.

Note: The medium requires inoculation of a large amount of feces than in case of other selective media.

Storage

Keep dry at room temperature.

Three (3) years after manufacturing.

Further Information

See Page 53

Blue Light Broth

Rapid test for Coliform group and *E. coli*

Code 1000211 – 300g (Granule)

Formula (in 17.4g for 1 liter)

Peptone.....	5.0g	Sodium Chloride.....	5.0g
Sodium Pyruvate.....	1.0g	Sodium Lauryl Sulfate.....	0.1g
Dipotassium Phosphate.....	4.0g	Monopotassium Phosphate.....	1.0g
Potassium nitrate.....	1.0g		
Isopropyl- β -D(-)-Thiogalactopyranoside (IPTG).....			0.1g
5-bromo-4-chloro-3-indolyl- β -D-galactopyranoside (X-GAL).....			0.1g
4-Methyl-Umbelliferyl- β -D-Glucuronide (MUG).....			0.1g
pH 7.1 \pm 0.2			

Directions

Suspend 17.4g of the dehydrated medium in 1,000 mL of distilled water; mix well and heat to dissolve the medium. Distribute necessary amounts of medium into test tube. Sterilize by autoclaving at 121°C for 15 minutes, and cool down rapidly for use.

Remarks

The medium is used for the test for Coliform group and *E. coli* in foodstuffs or water. Coliform group change the medium blue to blue-green after incubation at 35 – 37°C for 24 hours. Fluorescent is observed when the positive medium is exposed to UV light (365nm), if *E. coli* was present.

X-GAL (colorimetric enzyme substrate) in the medium is decomposed to bring out blue/blue-green color by β -galactosidase produced by Coliform group.

MUG in the medium is decomposed by β -Glucuronidase to be produced by *E. coli* to isolate a fluorescence substance of 4-Methyl-Umbelliferone.

Storage

Keep dry at room temperature.

Three (3) years after manufacturing.

Further Information

See Page 53

MacConkey Agar

MacConkey Agar, Modified 2

Code 05036– 300g

Formula (in 50.0g for 1 liter)

Peptone.....	19.0g	Crystal Violet.....	0.001g
Lactose.....	10.0g	Neutral Red.....	0.03g
Bile Salts.....	1.0g	Agar.....	15.0g
Sodium Chloride.....	5.0g		
pH 7.2±0.1			

Directions

Suspend 50.0g of the dehydrated medium in 1,000 ml of distilled water, mix well and heat to dissolve the medium. Sterilize by autoclaving at 121°C for 15 minutes, and distribute about 20ml into Petri dishes aseptically, and dry the surface of the medium before use.

Determinations

Incubate at 37°C for 18 – 20 hours.

Pathogenic bacteria form semitransparent colonies.

Colonies of *E. coli* assume an intensely redbrick color and produce light pink or redbrick sediments around.

Remarks

Features:

1. Differentiation between pathogenic bacteria and nonpathogenic bacteria on the medium is easier than in case of Endo's or Drigalski's medium.
2. Even some organisms that cannot grow on SS Agar grow on the medium.

Note:

Aseptic procedures should be taken for this medium contrary to those of SS Agar.

Storage

Keep dry at room temperature.

Three (3) years after manufacturing.

Further Information

See Page 53

MacConkey Sorbitol Agar

Isolation medium for *E. coli* O157

Code 05643 – 300g

Formula (in 50.0g for 1 liter)

Bile Salt.....	1.0g
Peptone.....	19.0g
D-Sorbitol	10.0g
Sodium Chloride.....	5.0g
Crystal Violet.....	0.001g
Neutral Red.....	0.03g
Agar.....	15.0g
pH 7.2±0.1	

Directions

Suspend 50.0g of the dehydrated medium in 1,000 mL of distilled water; mix well and heat to dissolve the medium. Sterilize by autoclaving at 121°C for 15 minutes. Distribute necessary amounts of medium into Petri dish.

Remarks

The medium is used for an isolation of *E. coli* O157 based on the nature that O157 does not ferment sorbitol while other *E. coli* do. *E. coli* O157 develops a semitransparent colony after incubation at 37°C for 18 – 20 hours. For the meanwhile other *E. coli* develop a redbrick color colony with light peach or redbrick zone of desoxycholic acid precipitated.

Growths of gram positive cocci are inhibited.

Storage

Keep dry at room temperature.
Three (3) years after manufacturing.

Further Information

See Page 53

Lauryl Sulfate MUG Broth

Rapid detection of Coliform group and *E. coli*

Code 05639 – 300g

Formula (in 35.7g for 1 liter)

Peptone.....	20.0g	Lactose.....	5.0g
Sodium Chloride.....	5.0g	Sodium Lauryl Sulfate.....	0.1g
Dipotassium Phosphate.....	2.75g	Monopotassium Phosphate.....	2.75g
4-Methyl-Umbelliferyl- β -D-Glucuronide (MUG).....	0.1g		
pH 6.8 ± 0.1			

Directions

Suspend 35.7g of the dehydrated medium in 1,000 mL of distilled water; mix well and heat to dissolve the medium. Distribute 10mL of medium into middle size of test tube that contains a fermentation tube or a Durham's tube. Sterilize by autoclaving at 121°C for 15 minutes, and cool down rapidly for use. Do not use any tubes that contain bubbles in a fermentation tube or a Durham's tube.

Remarks

The medium is used for the test for coliform group and *E. coli* in foodstuffs or water. Coliform group generate a gas after incubation at 35 – 37°C for 16 – 24 hours. Fluorescent is observed when the positive medium is exposed to UV light (365nm), if *E. coli* is presented. MUG in the medium is decomposed by β -Glucuronidase to be produced by *E. coli* to isolate a fluorescence substance of 4-Methyl-Umbelliferone.

Storage

Keep dry at room temperature.
Three (3) years after manufacturing.

Further Information

See Page 53

X-GAL Agar

Chromogenic Enzyme substrate agar for the test of Coliform group

Code 05642 – 300g

Formula (in 44.3g for 1 liter)

Peptone.....	15.0g	Yeast Extract.....	5.0g
Sodium Pyruvate.....	1.0g	Sodium Chloride.....	5.0g
Disodium Phosphate.....	2.0g	Potassium nitrate.....	1.0g
Sodium Lauryl Sulfate.....	0.15g	Agar.....	15.0g
5-bromo-4-chloro-3-indolyl- β -D-galactopyranoside (X-GAL).....	0.15g		
pH 7.1 \pm 0.2			

Directions

Suspend 44.3g of the dehydrated medium in 1,000 mL of distilled water; mix well and heat to dissolve the medium. Sterilize by autoclaving at 121°C for 15 minutes. Distribute necessary amounts of medium into Petri dish.

For the mixing and dilution culture of liquid specimen, the medium may be used in pour plate technique.

Remarks

The medium is used for the test for coliform group in foods and drinks. Coliform group develop a blue to blue-green colony after incubation at 35 – 37°C for 18 - 22 hours. Growths of almost all other bacteria are inhibited, and they develop only white colony even if they grow.

X-GAL (colorimetric enzyme substrate) in the medium is decomposed to bring out blue/blue-green color by β -galactosidase produced by coliform group.

Storage

Keep dry at room temperature.

Three (3) years after manufacturing.

Further Information

See Page 53

XM-G Agar

Chromogenic Enzyme substrate agar for the test of Coliform group and *E. coli*

Code 05647 – 300g

Formula (in 39.3g for 1 liter)

Peptone.....	10.0g	Sodium Pyruvate.....	1.0g
L-Tryptophan.....	1.0g	D-Sorbitol.....	1.0g
Sodium Chloride.....	5.0g	Monopotassium Phosphate.....	2.2g
Dipotassium Phosphate.....	2.7g	Potassium nitrate.....	1.0g
Sodium Lauryl Sulfate.....	0.2g	Agar.....	15.0g
5-bromo-4-chloro-3-indolyl- β -D-glucuronide (X-GLUC).....	0.1g		
5-bromo-6-chloro-3-indolyl- β -D-galactopyranoside (MAGENTA-GAL).....	0.1g		
pH 7.0 \pm 0.2			

Directions

Suspend 39.3g of the dehydrated medium in 1,000 mL of distilled water; mix well and heat to dissolve the medium. Sterilize by autoclaving at 121°C for 15 minutes. Distribute necessary amounts of medium into Petri dish.

Remarks

The medium is used for the screening test for *E. coli* and coliform group. *E. coli* has a β -glucuronidase that decomposes the substrate of X-GLUC to produce blue dye. For the meanwhile Coliform group have a β -galactosidase that decomposes the substrate of MAGENTA-GAL to produce red dye. *E. coli* has both enzymes, and may produce a blue to blue purple color.

Since *E. coli* O157 does not have a β -glucuronidase, they are identified as a Coliform group. Incubate at 35°C for 20 \pm 2 hours. Over-time incubation may foster growth of microorganism other than *E. coli* and Coliform group. Red color may be observed if the sample contains lactobacilli that has also a β -galactosidase.

Storage

Keep dry at room temperature.
Three (3) years after manufacturing.

Further Information

See Page 53

Staphylococcus Medium No. 110

Selective isolation medium for Staphylococci

Code 05234 – 300g

Formula (in 150.0g for 1 liter)

Peptone.....	10.0g	Sodium Chloride.....	75.0g
Yeast Extract.....	2.5g	Lactose.....	2.0g
Gelatin.....	30.0g	Dipotassium Phosphate.....	5.0g
Mannitol.....	10.0g	Agar.....	15.0g
pH 7.5±0.1			

]

Directions

Suspend 150.0g of the dehydrated medium in 1,000 ml of distilled water; mix well and heat to dissolve the medium. Sterilize by autoclaving at 121°C for 15 minutes. Distribute about 20ml of the medium into Petri dishes. Dry the surface of the plate before use.

Determinations

Incubate at 30°C for 48 hours or 37°C for 43 hours.

Pathogenic staphylococci form yellowish or lemon-colored colonies on the medium and have abilities to utilize mannitol, to liquefy gelatin (positive to Stone reaction), to coagulate plasma and to hemolyze rabbit blood.

Remarks

1. Most of bacteria except staphylococci are inhibited on the medium, and some of the non-pathogenic staphylococci are also inhibited to some extent.
2. Pathogenic staphylococci generally form yellowish or lemon-colored colonies on the medium, while non-pathogenic strains form white ones.
3. The ability for producing pigment, utilizing mannitol and liquefying gelatin can be directly examined on the medium.

Storage

Keep dry at room temperature.

Three (3) years after manufacturing.

Further Information

See Page 53

Mannitol Salt Agar

For isolation of *Staphylococcus*

Code 05236 –300g (granule)

Formula (in 111.0g for 1 liter)

Beef Extract.....	1.0g	D-Mannitol.....	10.0g
Peptone.....	10.0g	Phenol Red.....	0.025g
Sodium Chloride.....	75.0g	Agar.....	15.0g
pH 7.55±0.1			

Directions

Suspend 111.0g of the dehydrated medium in 1,000 ml of distilled water; mix well and heat to dissolve the medium. Sterilize by autoclaving at 121°C for 15 minutes. Distribute about 20ml of the medium into Petri dishes. Dry the surface of the plate before use.

Determinations

Incubate at 30°C for 48 hours or 37°C for 36 hours.

When staphylococci, which utilize mannitol, grow on the medium, they form yellow colonies and turn the medium around the colonies yellow, and *Staphylococcus aureus* can easily be differentiated from non-pathogenic staphylococci and other bacteria that do not utilize mannitol and form red colonies.

Remarks

As the medium contains salt of high concentration, non-halophytic bacteria are inhibited to grow, and staphylococci are selectively isolated. Checking the ability to utilize mannitol performs the screening test of pathogenic staphylococci

For differential medium of coagulase positive staphylococci, add egg yolk solution to the medium. To prepare the 10% salt egg yolk agar medium, suspend 111.0g of the dehydrated medium in 900 ml of distilled water, sterilize by autoclaving and cool down to 50°C, to which add 100 ml of egg yolk solution (Mixture of about 20g of egg yolk and 80 ml sterilized physiological saline solution).

Storage

Keep dry at room temperature.

Three (3) years after manufacturing.

Further Information

See Page 53

For Isolation of *Enterococcus***EF Agar Base**

(Appended reagent: TTC)

Code 05679 – 100g**Formula (in 54.5g for 1 liter)**

Calf Brain Extract Powder.....	8.5g	Beef Heart Extract Powder....	8.5g
Peptone.....	10.0g	Dextrose.....	10.0g
Dipotassium Phosphate.....	2.5g	Brom Thymol Blue.....	0.032g
Agar.....	15.0g		
2,3,5-Triphenyltetrazolium Chloride (TTC) (Appended reagent).....	0.15g		
pH 7.2±0.1			

Directions

Suspend 54.5g of the dehydrated medium together with 0.25g of Sodium azide in 1,000 mL of distilled water; mix well and heat to dissolve the medium. Keep at around 60°C, and then add 10mL of 1.5% TTC solution. Mix well and distribute the medium into Petri dishes.

1.5% TTC solution: Dissolve 0.15g of TTC (appended to the medium) in 10mL of distilled water, and store in brown bottle in dark and cool place.

Remarks

The medium is used for an isolation of *Enterococcus* in faces, water, milk, dairy products and processed meat.

Incubate at 35 – 37°C for 48 hours. *E. faecalis* forms pink or red brown colonies 0.5 – 2mm in diameter, while *E. faecium* forms yellow colonies. Growths of gram-negative bacilli are inhibited.

Since *Enterococcus* is an indicator of fecal pollution as in case of Coliform, the detection of it is important for the survey of fecal pollution.

Storage

Keep dry at room temperature. Three (3) years after manufacturing.

Further Information

See Page 53

NGKG Agar Base

For selective isolation of *Bacillus cereus*

Code 05282 – 300g (granule)

Formula (in 26.5g for 1 liter)

Peptone.....	1.0g	Polymyxin B Sulfate.....	50,000unit
Yeast Extract.....	0.5g	Phenol Red.....	0.025g
Sodium Chloride	4.0g	Agar.....	18.0g
Glycine.....	3.0g		
pH 6.8±0.1			

Directions

Suspend 26.5g of the dehydrated medium in 900ml of distilled water; mix well and heat to dissolve the medium. Sterilize by autoclaving at 121°C for 15 minutes. Maintain the medium at about 50°C, add 100ml of 20% egg-yolk suspension, mix well and distribute about 20ml amounts into Petri dishes.

For the preparation of egg-yolk suspension, add 20ml of egg-yolk to 80ml of sterilized saline solution aseptically, and mix well.

Determinations

After drying the surface of the plate, smear the specimens and incubate at 30°C for 18 – 24 hours.

B. cereus forms white and slightly thick colonies with an irregular margin on the medium and shows the lecithinase reaction. In case of the lecithinase reaction, their colonies form a zone of opacity, and the medium around the colonies presents a red color. The growth of miscellaneous bacteria except *B. cereus* may be inhibited. Even if some bacteria grow, their colonies are small and do not show the lecithinase reaction.

Remarks

The medium was designed for the selection and detection of *B. cereus* from the contaminated foods.

On the medium, the spores of *B. cereus* are formed well. After incubation at 30°C for 18 hours, the sporulation can be determined by observing them microscopically.

Storage

Keep dry at room temperature.

Three (3) years after manufacturing.

Further Information

See Page 53

GAM Agar, Modified

(For common anaerobic culture and susceptibility test)
(Gifu Anaerobic Medium Agar, Modified)

Code 05426 – 300g

Formula (in 56.7g for 1 liter)

Peptone.....	5.0g	L-Tryptophan.....	0.2g
Soya Peptone.....	3.0g	Monopotassium Phosphate.....	2.5g
Proteose Peptone.....	5.0g	Sodium Chloride.....	3.0g
Digested Serum.....	10.0g	L-Cysteine Hydrochloride.....	0.3g
Yeast Extract.....	2.5g	Sodium Thioglycollate.....	0.3g
Beef Extract.....	2.2g	L-Arginine.....	1.0g
Liver Extract.....	1.2g	Vitamin K ₁	5mg
Dextrose.....	0.5g	Hemin.....	10mg
Soluble Starch.....	5.0g	Agar.....	15.0g
pH 7.3±0.1			

Directions

Suspend 56.7g of the dehydrated medium in 1,000 ml of distilled water, and heat to dissolve the medium. Sterilize by autoclaving at 115°C for 15 minutes. Distribute the medium into Petri dishes.

Prepared medium should be used in a day of preparation, or keep in anaerobic condition.

Remark

The medium is a modified medium of GAM Agar that is developed by Medical School of Gifu University, Japan to isolate and culture anaerobic bacteria from clinical specimens. The medium is also used for susceptibility tests other than sulfa drugs.

Color of the medium is more pale and better transparency

Storage

Keep dry at room temperature.

Three (3) years after manufacturing.

Further Information

See Page 53

GAM Broth, Modified

(For common anaerobic culture and susceptibility test)
(Gifu Anaerobic Medium Broth, Modified)

Code 05433 –100g

Formula (in 41.7g for 1 liter)

Peptone.....	5.0g	L-Tryptophan.....	0.2g
Soya Peptone.....	3.0g	Monopotassium Phosphate.....	2.5g
Proteose Peptone.....	5.0g	Sodium Chloride.....	3.0g
Digested Serum.....	10.0g	L-Cysteine Hydrochloride.....	0.3g
Yeast Extract.....	2.5g	Sodium Thioglycollate.....	0.3g
Beef Extract.....	2.2g	L-Arginine.....	1.0g
Liver Extract.....	1.2g	Vitamin K ₁	5mg
Dextrose.....	0.5g	Hemin.....	10mg
Soluble Starch.....	5.0g		
pH 7.3±0.1			

Directions

Suspend 41.7g of the dehydrated medium in 1,000 ml of distilled water, and heat to dissolve the medium. Distribute the medium into appropriate containers. Sterilize by autoclaving at 115°C for 15 minutes, and cool down quickly (Do not shake the medium!).

Remark

The composition of the medium is exactly same as GAM Agar, Modified excluding agar. The medium is liquid and used for isolation and cultivation of anaerobic bacteria from clinical specimens. The medium is also used for susceptibility tests other than sulfa drugs, in particular for the micro liquid dilution method of anaerobic bacteria.

Storage

Keep dry at room temperature.
Three (3) years after manufacturing.

Further Information

See Page 53

Bacteroides Agar

For differentiation and selective isolation of *Fusobacterium* and *Bacteroides*

Code 05440 – 100g

Formula (in 74.0g for 1 liter)

Peptone.....	20.0g	Monopotassium Phosphate.....	2.5g
Soya Peptone.....	1.5g	Sodium Chloride.....	3.0g
Digested Serum Powder.....	6.75g	Soluble Starch.....	5.0g
Liver Extract.....	0.6g	L-Cysteine Hydrochloride.....	0.3g
Meat Extract.....	6.15g	Sodium Thioglycollate.....	0.3g
Yeast Extract.....	10.0g	Colistin.....	1,000,000 units
Hemin.....	0.003g	Neomycin.....	0.2g
Dextrose.....	3.0g	Brilliant Green.....	0.001g
Agar.....	14.7g	pH 7.1 ± 0.1	

Directions

Suspend 74.0g of the dehydrated medium in 1,000ml of distilled water; mix well and heat to dissolve the medium. Sterilize by autoclaving at 121°C for 15 minutes, and distribute the medium into Petri dishes. Dry the surface sufficiently before use. It is desirable to use the plate within 3 – 5 hours after preparation.

Determinations

Differentiation: The medium is able to identify *Bacteroides* from the organisms that are isolated from clinical specimens and confirmed to be asporogenic gram-negative bacilli of obligate anaerobes. Inoculate the bacilli into GAM Semisolid and incubate for 24 – 48 hours. Take the organisms, streak the surface of the plate, and incubate at 37°C for 24 – 48 hours under an anaerobic condition. The organisms are identified as *Bacteroides*, when they grow well on this plate.

Selective isolation: When clinical specimens are considered to be lightly contaminated with other organisms, *Bacteroides* can be isolated selectively by smearing them directly. (*Fusobacterium* does not grow on this medium, but some aerobic cocci may grow on it). For the quantitative culture, dilute a sample with the following solution (prepared by 1/15M phosphate-buffered solution, pH 7.2, containing 1g of polysorbate 80, 1g of L-cysteine hydrochloride and 1g of agar per liter by autoclaving at 110°C for 15 minutes).

Remarks

The medium was devised for the differentiation of *Bacteroides*. Colistin, neomycin and brilliant green have no influence on the growth of *Bacteroides* but inhibit *Fusobacterium* and other bacteria. Asporogenic gram-negative bacilli taken from the clinical specimens are mainly *Bacteroides* or *Fusobacterium*. The medium also permits the selective growth of *Bacteroides* from the clinical specimens.

Storage

Keep dry at room temperature.
Three (3) years after manufacturing.

Further Information

See Page 53

FM Agar, Modified

For differentiation and selective isolation of *Fusobacterium*

Code 05441 – 100g

Formula (in 74.0g for 1 liter)

Peptone.....	20.0g	Monopotassium Phosphate.....	2.5g
Soya Peptone.....	1.5g	Sodium Chloride.....	3.0g
Digested Serum Powder.....	6.75g	Soluble Starch.....	5.0g
Liver Extract.....	0.6g	L-Cysteine Hydrochloride.....	0.3g
Meat Extract.....	6.15g	Sodium Thioglycollate.....	0.3g
Yeast Extract.....	10.0g	Neomycin.....	0.2g
Dextrose.....	3.0g	Crystal Violet.....	0.01g
Agar.....	14.7g		
pH 7.1±0.1			

Directions

Suspend 74.0g of the dehydrated medium in 1,000ml of distilled water, and heat to dissolve the medium. Do not sterilize by autoclaving. Distribute the medium into Petri dishes. It is recommended to use the medium in 3 – 5 hours after preparation.

Determination

The medium is designed to differentiate *Fusobacterium*. Neomycin and crystal violet have no influence on the growth of *Fusobacterium*, but inhibit *Bacteroides* and other bacteria.

Since it is well known fact that almost all gram-negative rods, which do not form spore, are *Fusobacterium* and *Bacteroides*. the medium is used for the isolation of *Fusobacterium*. Among of non-spore forming, gram-negative anaerobic rods, and the bacteria grown on the medium is to be identified as *Fusobacterium*.

Remarks

To prepare the medium, heating should be minimum, and use the medium in 3 – 5 hours after preparation.

Incubate the medium avoiding excess disclosure to the air.

Storage

Keep dry at room temperature.

Three (3) years after manufacturing.

Further Information

See Page 53

ANAERO MATE-P ANAERO MATE-J

Simple Anaerobic Culture System

**Code Anaero Mate-P – 06593 (For 100 plates: 50 set x for 2 plates)
 Anaero Mate-J – 06594 (For 500 plates: 50 set x for 10 plates)**

Set component

Anaerobic Culture Bag.....50 pcs
Anaerobic Culture Agent (absorb oxygen and generate carbon dioxide).....50 pcs
Oxygen Detector50 pcs
Seal Bar.....2 pcs

Directions

Set petri dishes (plates) in Anaerobic Culture Bag, and then put Anaerobic Culture Agent and Oxygen detector in the bag subsequently. Seal the bag with Seal Bar immediately and put the whole bag(s) in incubator for cultivation.

Storage

Keep at 2 – 8°C. Do not freeze.
One year (12 months) after manufacturing.

Further Information

See Page 53

Heart Infusion Agar

For isolation of general bacteria

Code 05503– 300g (Granule)

Formula (in 40.0g for 1 liter)

Heart Extract Powder.....	10.0g	Sodium Chloride.....	5.0g
Peptone.....	10.0g	Agar.....	15.0g
pH 7.4 ± 0.1			

Directions

Suspend 40.0g of the dehydrated medium in 1,000 ml of distilled water; mix well and heat to dissolve the medium. Sterilize by autoclaving at 121°C for 15 minutes. Distribute into Petri dishes or test tubes according to the purpose.

Remarks

Since the medium has an excellent ability to support the growth of bacteria, it is suitable for the cultivation of many fastidious pathogenic bacteria.

Since *E. coli*, *Salmonella* and *Shigella*, which frequently mutate, are stable during the preservation on this medium, it is useful for the test of their antigenicity and virulence.

Storage

Keep dry at room temperature.

Three (3) years after manufacturing.

Further Information

See Page 53

Heart Infusion Broth

For isolation of general bacteria

Code 05505– 100g

Formula (in 25.0g for 1 liter)

Heart Extract Powder.....	10.0g	Sodium Chloride.....	5.0g
Peptone.....	10.0g		
pH	7.2 ± 0.1		

Directions

Suspend 25.0g of the dehydrated medium in 1,000 ml of distilled water; mix well and heat to dissolve the medium. Distribute into adequate containers according to the purpose, and sterilize by autoclaving at 121°C for 15 minutes.

Remarks

Since the medium has an excellent ability to support the growth of many kinds of bacteria, it is suitable for the cultivation and preservation of many fastidious bacteria.

The medium does not affect the biological nature of bacteria, and it is also suitable for precise biological tests.

The medium may be employed in the routine cultivation and for various other purposes, because the medium does not contain any carbohydrate such as dextrose.

Storage

Keep dry at room temperature.

Three (3) years after manufacturing.

Further Information

See Page 53

Brain Heart Infusion Agar

For isolation of general bacteria

Code 05506– 300g

Formula (in 50.0g for 1 liter)

Calf Brain Extract Powder.....	7.5g	Heart Extract Powder.....	8.0g
Peptone.....	10.0g	Dextrose.....	2.0g
Sodium Chloride.....	5.0g	Dipotassium Phosphate.....	2.5g
Agar.....	15.0g		
pH	7.2±0.1		

Directions

Suspend 50.0g of the dehydrated medium in 1,000 ml of distilled water; mix well and heat to dissolve the medium. Sterilize by autoclaving at 121°C for 15 minutes. Use as plates or slants.

Remarks

Since the medium has an excellent ability to support the growth of many kinds of bacteria, it is suitable for the cultivation and preservation of many nutritionally fastidious bacteria.

The medium does not affect the biological nature of bacteria, and it is suitable for precise biological tests.

Storage

Keep dry at room temperature.
Three (3) years after manufacturing.

Further Information

See Page 53

Brain Heart Infusion Broth

For isolation of general bacteria

Code 05508– 300g

Formula (in 35.0g for 1 liter)

Calf Brain Extract Powder.....	7.5g	Heart Extract Powder.....	8.0g
Peptone.....	10.0g	Dextrose.....	2.0g
Sodium Chloride.....	5.0g	Dipotassium Phosphate.....	2.5g
pH 7.2±0.1			

Directions

Suspend 35.0g of the dehydrated medium in 1,000 ml of distilled water; mix well and heat to dissolve the medium. Sterilize by autoclaving at 121°C for 15 minutes.

Remarks

Since the medium has an excellent ability to support the growth of many kinds of bacteria, it is suitable for the cultivation and preservation of many nutritionally fastidious bacteria. The medium is also used for the cultivation of blood specimens.

The medium does not affect the biological nature of bacteria, and it is suitable for precise biological tests.

Storage

Keep dry at room temperature.
Three (3) years after manufacturing.

Further Information

See Page 53

Trypto-Soya Agar

SCD Agar

For isolation of general bacteria

Code 05516 – 300g (Granule)

Formula (in 40.0g for 1 liter)

Peptone.....	15.0g	Sodium Chloride.....	5.0g
Soya Peptone.....	5.0g	Agar.....	15.0g
pH 7.2±0.1			

Directions

Suspend 40.0g of the dehydrated medium in 1,000 ml of distilled water; mix well and heat to dissolve the medium. Sterilize by autoclaving at 121°C for 15 minutes. Use as plates or slants.

Remarks

The medium is fit for various purposes because it permits the growth of fastidious bacteria that do grow on Nutrient Agar or Heart Infusion Agar, and is markedly superior to conventional media as a basal medium of blood agar.

Storage

Keep dry at room temperature.
Three (3) years after manufacturing.

Further Information

See Page 53

Mueller-Hinton Agar-N

For isolation of general bacteria

Code 05533 – 300g

Formula (in 38.0g for 1 liter)

Beef extract	2.0g	Casamino acid	17.5g
Soluble starch.....	1.5g	Agar.....	17.0g
pH 7.3±0.1			

Directions

Suspend 38.0g of the dehydrated medium in 1,000 ml of distilled water; mix well and heat to dissolve the medium. Sterilize by autoclaving at 121°C for 15 minutes. Use for plates and thickness of medium shall be 4mm (for 9cm plate, pour 25mL of autoclaved medium).

Remarks

This medium is also fit for Anti-susceptibility testing which applicable for K-B method (following to NCCLS document M2-A6).

Storage

Keep dry at room temperature.

Three (3) years after manufacturing.

Further Information

See Page 53

Sabouraud Agar

For isolation of Yeast and Mold

Code 05701 – 300g

Formula (in 65.0g for 1 liter)

Peptone.....	10.0g	Agar.....	15.0g
Dextrose.....	40.0g		
pH 5.9 ± 0.2			

Directions

Suspend 65.0g of the dehydrated medium in 1,000 mL of distilled water; mix well and heat to dissolve the medium. Sterilize by autoclaving at 121°C for 15 minutes, and use as a plate or a slant.

Determinations

Incubate at 37°C or at 25°C. For the isolation of *Candida*, incubate for 48 – 72 hours and in case of other fungi, incubation for 5 days may sometimes be necessary. *Candida* forms colonies with white, opaque, and wet swelling, which gradually change to light brown.

Remarks

Growth of *Candida* on the medium is so well that it can survive for a long period without morphological variations. As the medium cannot inhibit the contaminant microbes, it is desirable to use the selective media such as Mycobiotic Agar and Candida GE agar simultaneously for the isolation.

Storage

Keep dry at room temperature.
Three (3) years after manufacturing.

Further Information

See Page 53

Candida GE Agar

For isolation of *Candida*

Code 05703 – 100g

Formula (in 62.0g for 1 liter)

Yeast Extract.....	10.0g	Nitrofurantoin.....	0.5g
Peptone.....	8.5g	Agar.....	13.0g
Dextrose.....	30.0g		
pH 6.0 ± 0.1			

Directions

Suspend 62.0g of the dehydrated medium in 1,000mL of distilled water; mix well and heat to dissolve the medium. Distribute about 15 – 20mL amounts into Petri dishes. No sterilizing procedure such as autoclaving is needed.

Determinations

After drying the surface, inoculate the slime in oral or vaginal secretions with a sterilized swab and spread on the surface of the medium. Incubate at 35 – 37°C for 2 – 3 days.

Candida forms specific colonies that are round or oval in shape with a diameter of 3 – 5mm and are wet and opaque, and have specific odor and luster.

C. albicans forms colonies of creamy or off-white color after 48 – 72 hours incubation, which gradually turn pale brownish.

C. krusei forms flat and irregularly shaped colonies with no luster, while the other *Candida* form brown or pale brown colonies.

Remarks

Features:

1. The medium inhibits the growth of miscellaneous bacteria and permits the growth of *Candida* specifically.
2. Nitrofurantoin markedly or completely inhibits gram-negative bacilli. Therefore the medium facilitates the detection of *Candida*.

Note:

It is necessary to smear the specimens heavily on the medium.

Storage

Keep dry at room temperature.

Three (3) years after manufacturing.

Further Information

See Page 53

ID Test EB-20

Identification Test Kit for Glucose fermentative Gram Negative bacilli

Code 06626 – 25 tests

06628 – ID Test EB-20 EB Reagent 100 tests

06615 – ID Test EB-20 EB Plate 100 tests (additional plates)

06616 – ID Test EB-20 EB Broth 100 tests (additional broth)

Kit component (06626)

ID Test EB Plate.....25 plates x 20 wells
Sterilized Liquid Paraffin.....20mL x 1 vial
Result Form.....13 sheets x 2 sets

ID Test EB Broth.....2.5mL x 25 tubes
Color Chart.....1 sheet

Kit component (06628) for 100 tests

10% Ferric Chloride Test solution..... 6mL x 1 vial
Kovacs' Reagent..... 6mL x 1 vial
6% α -naphthol Test solution.....6mL x 1 vial
40% Potassium Hydroxide Test solution... 6mL x 1 vial

Features of ID Test EB-20

Identification of bacteria is very complicated and there are so many biochemical test need to be done. ID Test EB-20 is based on the theory of the numerical classification method, which could receive the combination with 20 kinds of biochemical reactions.

Only small amount of sample (0.1mL) is needed for each test item. Also the kit is stable for 1 year at room temperature.

Purpose of ID Test EB-20

Identification of Glucose fermentative Gram Negative bacilli

Directions for operation of the kit

1. Streak samples onto blood agar or other isolation medium. After incubation of the medium, confirm colonies will be tested are Gram Negative bacilli and Glucose fermentative.
2. The test target bacteria shall be proliferated on enrichment medium.
3. Fish colonies from enrichment medium for adjustment the concentration (suspension) that corresponds to No.1 McFarland turbidity standard.
4. Inoculate 0.1mL of above suspension into each well (20 wells). Layer 3-5 drops of sterilized liquid paraffin on 8 certain wells and incubate (place the cover on the plate and incubate at 37°C for 18 – 20 hours).
5. After incubation, add 1 drop of 10% Ferric Chloride Test solution (PPA), Kovacs' Reagent (IND), 6% α -naphthol Test solution and 40% Potassium Hydroxide Test solution (VP) into certain wells.
6. After addition of specific reagents, compare the color developed with the attached color chart and determine as positive or negative results. Simultaneously, an oxidase test shall be performed and record in the result form.

Interpretations

Refer ID Test EB-20 ANALYTICAL PROFILE for identification.

Storage and shelf life

Code 06626: Keep at Room Temperature. One year (12 months) after manufacturing.

Code 06628: Keep at 2 – 10°C. One year (12 months) after manufacturing.

Further Information

See Page 53

ID Test NF-18

Identification Test Kit for Glucose Non-fermentative Gram Negative bacilli

Code 06629 – 25 tests

06631 – ID Test NF-18 NF Reagent 100 tests

06617 – ID Test NF-18 NF Plate 100 tests (additional plates)

06618 – ID Test NF-18 NF Broth 100 tests (additional broth)

Kit component (06629)

ID Test NF Plate.....25 plates x 18 wells
Sterilized Liquid Paraffin.....20mL x 1 vial
Result Form.....13 sheets x 2 sets

ID Test NF Broth.....2.5mL x 25 tubes
Color Chart.....1 sheet

Kit component (06631) for 100 tests

Kovacs' Reagent.....6mL x 1 vial
Nitrite Detection Test Solution I.....6mL x 1 vial
Nitrite Detection Test Solution II.....6mL x 1 vial
Zinc Powder.....6mL x 1 vial

Features of ID Test NF-18

Identification of bacteria is very complicated and there are so many biochemical test need to be done. ID Test NF-18 is based on the theory of the numerical classification method, which could receive the combination with 18 kinds of biochemical reactions.

Only small amount of sample (0.1mL) is needed for each test item. Also the kit is stable for 1 year at room temperature.

Purpose of ID Test NF-18

Identification of Glucose Non-fermentative Gram Negative bacilli

Directions for operation of the kit

1. Streak samples onto blood agar or other isolation medium. After incubation of the medium, confirm colonies will be tested are Gram Negative bacilli and Glucose fermentative.
2. The test target bacteria shall be proliferated on enrichment medium.
3. Fish colonies from enrichment medium for adjustment the concentration (suspension) that corresponds to No.0.5 McFarland turbidity standard.
4. Inoculate 0.1mL of above suspension into each well (18 wells). Layer 3-5 drops of sterilized liquid paraffin on 9 certain wells and incubate (place the cover on the plate and incubate at 30 or 37°C for 22 – 24 hours).
5. After incubation, add 1 drop of Kovacs' Reagent (IND), Nitrate Detection Test Solution (I and II) (NIT) into certain wells.
6. After addition of specific reagents, compare the color developed with the attached color chart and determine as positive or negative results. Simultaneously, an oxidase test shall be performed and record in the result form.

Interpretations

Refer ID Test NF-18 ANALYTICAL PROFILE for identification.

Storage and shelf life

Code 06629: Keep at Room Temperature. One year (12 months) after manufacturing.

Code 06631: Keep at 2 – 10°C. One year (12 months) after manufacturing.

Further Information

See Page 53

N-ID Test SP-18

Identification Test Kit for Staphylococci

Code 06637 – 25 tests
06638 – N-ID Test SP-18 SP Reagent 100 tests
06613 – N-ID Test SP-18 SP Plate 100 tests (additional plates)
06614 – N-ID Test SP-18 SP Broth 100 tests (additional broth)

Kit component (06637)

N-ID Test SP Plate.....	25 plates x 18 wells	N-ID Test SP Broth.....	2.5mL x 25 tubes
Sterilized Liquid Paraffin.....	20mL x 1 vial	Color Chart.....	1 sheet
Result Form.....	13 sheets x 2 sets		

Kit component (06638) for 100 tests

6% α -naphthol Test solution.....	6mL x 1 vial
40% Potassium Hydroxide Test solution...	6mL x 1 vial
Nitrite Detection Test Solution I.....	6mL x 1 vial
Nitrite Detection Test Solution II.....	6mL x 1 vial

Features of N-ID Test SP-18

Identification of bacteria is very complicated and there are so many biochemical test need to be done. N-ID Test SP-18 is based on the theory of the numerical classification method, which could receive the combination with 18 kinds of biochemical reactions. Only small amount of sample (0.1mL) is needed for each test item. Also the kit is stable for 1 year at room temperature.

Purpose of ID Test SP-18

For identification of various species of Staphylococci.

Directions for operation of the kit

1. Streak samples onto blood agar or other isolation medium. After incubation of the medium, confirm colonies will be tested are Gram Negative bacilli and Glucose fermentative.
2. The test target bacteria shall be proliferated on enrichment medium.
3. Fish colonies from enrichment medium for adjustment the concentration (suspension) that corresponds to No.3 McFarland turbidity standard.
4. Inoculate 0.1mL of above suspension into each well (18 wells). Layer 3-5 drops of sterilized liquid paraffin on 2 certain wells and incubate (place the cover on the plate and incubate at 30 or 37°C for 22 – 24 hours. (Wells 1 and 2 are empty wells and inoculation of bacterial sample solution is not necessary.
5. After incubation, add 1 drop of Kovacs' Reagent (IND), 6% α -naphthol and 40% Potassium Hydroxide Test solution (VP) into certain wells.
6. After addition of specific reagents, compare the color developed with the attached color chart and determine as positive or negative results.

Interpretations

Refer N-ID Test SP-18 ANALYTICAL PROFILE for identification.

Storage and shelf life

Code 06637: Keep at Room Temperature. One year (12 months) after manufacturing.
 Code 06638: Keep at 2 – 10°C. One year (12 months) after manufacturing.

Further Information

Germany / Austria / Switzerland

HyServe GmbH & Co. KG Hechenrainer Strasse 24 82449 Uffing Germany Phone: 0049-8846-1344 Fax: 0049-8846-1342 info@hyserve.com www.HyServe.com	QMT Ingenieurbüro Poststraße 29 30890 Barsinghausen Deutschland Phone: 05105-5209841 Fax: 05105-5209841 info@qmti.de www.qmti.de
---	--